EnoGeneFec™ 2200 Transfection Reagent

EGF2200

EnoGeneFec™ 2200 Transfection Reagent

Catalog Number: EGF2200-500, EGF2200-1000, EGF2200-1500

Amount: 500µl 1000µl 1500µl

Storage/Stability: Store at 4°C/1 year

Introduction

EnoGeneFec™ 2200 Transfection Reagent comprises of unique formulations of polycations and liposomes representing a new class of transfection reagent designed for primary mammalian cells to get outstanding transfection results, which will guarantee higher transfection efficiency and lower cytotoxicity for almost all primary cell lines. EnoGeneFec™ 2200 Transfection Reagent offers significant advantages over many other transfection methods.

EnoGeneFec™ 2200 Transfection Reagent offers:

- High transfection efficiency of plasmids, antisense oligonucleotides, siRNA in almost all primary cell types
- Fast procedure — transfection-complex formation in just 20 minutes
- Decreased cytotoxicity

Kit Contents

Cat:	EnoGeneFec™ 2200 Transfection Reagent
EGF2200-500	500µL
EGF2200-1000	1.0mL
EGF2200-1500	2.0mL

Important Notes

Required amount of EnoGeneFec™ 2200 Transfection Reagent

For initial optimization experiments, transfect a monolayer of cells that is 70–90% confluent, using 2:1~ 5:1 ratios of EnoGeneFec™ 2200 Transfection Reagent (µl) to DNA (µg), respectively. For most cell types, the EnoGeneFec™ 2200 Transfection Reagent: DNA (2:1) ratio provides excellent transfection level. For some hard-to-transfect cells, the EnoGeneFec™ 2200 Transfection Reagent: DNA (5:1) ratio is recommended.

Cell culture

A healthy cell culture lays the foundation for successful transfection. Different cells or cell lines have very specific media, serum, and supplement requirements. Low passage number (<50 splitting cycles) ensures that the cell genotype has not become altered. We recommend subculturing cells 24 hours before transfection. This provides normal cell metabolism and increases the likelihood of DNA uptake.

Effect of serum

We do not recommend using serum during complex formation between EnoGeneFec™ 2200 Transfection Reagent and plasmid DNA, as serum may inhibit complex formation.

Plasmid DNA quality

It is critical to accurately determine the plasmid DNA concentration using 260nm absorption. DNA content must be determined by 260nm absorption (estimates of DNA content based on the intensity of gel bands are not sufficiently accurate). Determine the DNA purity using a 260 nm/280 nm ratio; the ratio should be 1.8.

Procedures

Use the following conditions as guidelines to transfect primary mammalian cells in a 6-well or 35mm dish format. For other culture vessels, please refer to Table 1.

Table 1: Reagent Quantities for Different Culture Vessels

Culture Vessel	Surface area per well (cm ²)	Typical total volume of medium (ml)	Starting mass of DNA (µg) in medium volume (µl)	EnoGeneFec™ 2200 in medium volume	Transfection medium vol.
96 well	0.3	0.1-0.2	0.2 µg in 20 µl	0.4-1.0 µl in 20 µl	50 µl
24 well	2	0.5	(0.2-)1.0 µg in 25 µl	2-5 µl in 25 µl	0.4 ml
12 well	4	1.0	(0.5-)2.0 µg in 100 µl	4-10 µl in 100 µl	0.6 ml
6 well	10	2.0	(1.0-)5.0 µg in 100 µl	10-25 µl in 100 µl	0.8 ml
35 mm	10	2.0	(1.0-)5.0 µg in 100 µl	10-25 µl in 100 µl	0.8 ml
60 mm	20	5.0	(3.0-)10.0 µg in 500 µl	20-50 µl in 500 µl	2.4 ml
10 cm	60	10.0	(8.0-)20.0 µg in 1500 µl	40-100 µl in 800µl	6.4 ml

For *in vitro* primary adherent cells

1. 18 to 24 hours prior to transfection, trypsinize and adjust the cell concentration, seed cells at a density of $1-3 \times 10^5$ cells per well in 2.0ml of appropriate growth medium (with serum and antibiotics if cells are cultured in the presence of them). Incubate the cells at 37°C in a CO₂ incubator until cells are 70% to 90% confluent at the time of transfection.

2. For each transfection sample, prepare the complexes as follows:

Solution A: Dilute 2.0 µg of DNA into 100 µl of serum-free, antibiotic-free medium.

Solution B: Vortex EnoGeneFec™ 2200 Transfection Reagent thoroughly prior use, then dilute 10 µl of EnoGeneFec™ 2200 Transfection Reagent in 100 µl serum-free, antibiotic-free medium.

Incubate Solution A and B at room temperature for 5 minutes.

3. Combine the solutions, mix gently to ensure uniform distribution and incubate for 20 minutes at room temperature. NOTE: *Complexes are stable at room temperature for 3-5 hours.*

4. Add 0.8 ml of serum-free, antibiotic-free medium to EnoGeneFec™-DNA complex. Mix solution gently.

5. Remove growth medium from the cells and add 1.0 ml of EnoGeneFec™-DNA solution to the each well containing cells.

6. After 5-8 hours, remove transfection solution and add 2.0 ml of the appropriate growth medium (with serum and antibiotics). Incubate the cells at 37°C in a CO₂ incubator for a total of 24-72 hours.

7. Assay cell extracts or phenotype for gene activity 24-72 hours after the start of transfection depending on the cell types and promoter activity. To make stable cell lines: Passage cells at a 1:10 (or higher dilution) into fresh growth medium 24 hours post transfection. Selection medium can be added the following day if desired.

For *in vitro* primary suspension cells

1. Just prior to preparing complexes, plate $3-5 \times 10^5$ suspension cells in 0.8ml of serum free medium without antibiotics. Since transfection efficiency is sensitive to culture confluence, it is important to maintain a standard seeding protocol from experiment to experiment.

2. For each transfection sample, prepare the complexes as follows:

Solution A: Dilute 2.0 µg of DNA into 100 µl of serum-free, antibiotic-free medium.

Solution B: Vortex EnoGeneFec™ 2100Transfection Reagent thoroughly prior use, then dilute 10 µl of EnoGeneFec™ 2100Transfection Reagent in 100 µl serum-free, antibiotic-free medium.

Incubate Solution A and B at room temperature for 5 minutes.

3. Combine the solutions, mix gently to ensure uniform distribution and incubate for 20 minutes at room temperature. NOTE: *Complexes are stable at room temperature for 3-5 hours.*

4. Add 0.2ml of the EnoGeneFec™-DNA solution into each well containing suspension cells in 0.8ml serum-free, antibiotic-free medium.
5. After 5-8 hours, add 0.1ml of FBS directly into each vessel.. Incubate the cells at 37°C in a CO₂ incubator for a total of 24-72 hours.
6. Assay cell extracts or phenotype for gene activity 24-72 hours after the start of transfection depending on the cell types and promoter activity. To make stable cell lines: Passage cells at a 1:10 (or higher dilution) into fresh growth medium 24 hours post transfection. Selection medium can be added the following day if desired.

Storage and stability

EnoGeneFec™ 2200 Transfection Reagent is provided in 1µg/µL concentration. It is shipped at room temperature and is stabilized for extended storage at +4°C for one year when very tightly closed.

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