EnoGeneFec™ 2300 Transfection Reagent

EGF2300

EnoGeneFec™ 2300 Transfection Reagent

Catalog Number: EGF2300-500, EGF2300-1000, EGF2300-1500

Amount: 500µl 1000µl 1500µl

Storage/Stability: Store at 4°C/1 year

Introduction

EnoGeneFec™ 2300 Transfection Reagent comprises of unique formulations of polycations and liposomes representing a new class of transfection reagent designed for outstanding *in vivo* transfection results, which will guarantee higher transfection efficiency and lower cytotoxicity. EnoGeneFec™ 2300 Transfection Reagent offers significant advantages over many other transfection methods.

EnoGeneFec™ 2300 Transfection Reagent offers:

- High transfection efficiency of plasmids, antisense oligonucleotides, siRNA
- Fast procedure — transfection-complex formation in just 20 minutes
- Decreased cytotoxicity and biodegradable
- Suitable for *in vivo* transfection

Kit Contents

| Cat: | EnoGeneFec™ 2300 Transfection Reagent |
|--------------|---------------------------------------|
| EGF2300-500 | 500µL |
| EGF2300-1000 | 1.0mL |
| EGF2300-1500 | 2.0mL |

Important Notes

Required amount of EnoGeneFec™ 2300 Transfection Reagent

For most cell types, the EnoGeneFec™ 2300 Transfection Reagent: DNA (2:1) ratio provides excellent transfection level. For some hard-to-transfect cells, the EnoGeneFec™ 2300 Transfection Reagent: DNA (5:1) ratio is recommended.

Effect of serum

We do not recommend using serum during complex formation between EnoGeneFec™ 2300 Transfection Reagent and plasmid DNA, as serum may inhibit complex formation.

Plasmid DNA quality

It is critical to accurately determine the plasmid DNA concentration using 260nm absorption. DNA content must be determined by 260nm absorption (estimates of DNA content based on the intensity of gel bands are not sufficiently accurate). Determine the DNA purity using a 260 nm/280 nm ratio; the ratio should be 1.8.

In vivo transfection

For *in vivo* transfection, disease animal model such as nude mice human cancer xenograft model shall be chose according to customer's research purpose. The *in vivo* dosage of oligos (including plasmids, antisense oligonucleotides, siRNA, etc) may be adjusted to the effective dosage based on customer's previous studies. For most condition, the EnoGeneFec™ 2300 Transfection Reagent (2:1 or 5:1) ratio provides excellent transfection level. The route of administration of EnoGeneFec™ 2300 Transfection Reagent: DNA complex may be designed according to future clinical application route such as *i.v.*, *i.m.*, *i.p.*, or intratumor injection, etc.

Procedures

Use the following conditions as guidelines for *in vivo* transfection.

The EnoGeneFec™ 2300 Transfection Reagent: DNA ratio is recommended according to the types of oligos (including plasmids, antisense oligonucleotides, siRNA, etc) as in Table 1 for customer's optimization.

Table 1: Reagent Quantities for Different Oligos

| Types of oligos | EnoGeneFec™ 2300 Transfection Reagent: DNA ratio |
|----------------------------|--|
| Plasmids | 2-5: 1 |
| Antisense oligonucleotides | 2-5: 1 |
| siRNA | 2-5: 1 |

For nude mice human cancer xenograft model gene therapy by intratumor injection as example:

1) Nude mice human cancer xenograft model preparation

Under SPF conditions, the tumor particles or Log phase cancer cells are inoculated s.c. of nude mice right axilla. The xenografted tumor was measured the length and width by vernier caliper. Nude mice are randomly divided into several groups when the bearing tumor reached 100 mm³ volume. Tumor length and width were measured twice every week to dynamic evaluate the anticancer activity of the tested compounds.

Tumor volume (TV) was calculated by the following formula:

$$TV = 1/2 \times a \times b^2$$

a and b represents length and width, respectively.

2) Preparation of EnoGeneFec™ 2300 Transfection Reagent: DNA complex

Solution A: Dilute 10.0 µg of DNA into 50 µl of serum-free, antibiotic-free medium.

Solution B: Vortex EnoGeneFec™ 2300 Transfection Reagent thoroughly prior use, then dilute 50 µl of EnoGeneFec™ 2300 Transfection Reagent in 50 µl serum-free, antibiotic-free medium.

Incubate Solution A and B at room temperature for 5 minutes.

Combine the solutions, mix gently to ensure uniform distribution and incubate for 20 minutes at room temperature. NOTE: *Complexes are stable at room temperature for 3-5 hours.*

3) Complex injection

The above 100 µl EnoGeneFec™ 2300 Transfection Reagent-DNA complex is injected directly into the tumor after degermination of the local tumor.

4) Multiple dosing

If multiple dosing is required, repeat step 2) and 3). After period of oligos DNA administration, nude mice are sacrificed and tumors are isolated by surgery for analysis.

Storage and stability

EnoGeneFec™ 2300 Transfection Reagent is provided in 1 µg/µL concentration. It is shipped at room temperature and is stabilized for extended storage at +4°C for one year when very tightly closed.

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